

Digital Slide Scanning Service

Send us your glass microscope slides and we'll scan them to create digital slides that can be viewed, annotated, shared, managed, and analyzed with ease. Our Digital Slide Scanning Service provides a simple, affordable, and fast way to convert glass slides into digital slides. Digitizing your glass slides is the first step toward digital pathology, a digital environment for the management and interpretation of digital slide information.

Let us convert your glass microscope slides into high-resolution, Internet-accessible digital slides.

Enable easy access to your digital slides by making them accessible from any workstation on the network, eliminating the delay of physically transporting glass slides from one location to another.

Leverage remote viewing capabilities as well as manage volumes of digital slides and supplemental data by adding Spectrum™ Plus Hosted, the most comprehensive, web-based digital pathology information management software. With Spectrum™ Plus, you will have a powerful application to make your digital slide information available anywhere, anytime over the Internet via secure web-based access. You'll benefit from an affordable and easy-to-manage digital slide hosting and management interface while we maintain the software, hardware, and other IT components for you.

Why wait when you can start taking advantage of digital pathology today?



KEY BENEFITS

- Minimize the task of managing glass slide sets
- Interactive whole-slide navigation superior to viewing static images
- Supplement digital slides with descriptive information by annotating regions of interest
- Distribute annotated slides with ease on DVD

HIGH-QUALITY DIGITAL SLIDES

- Entire slide captured at 20x or 40x objective
- Oil immersion scanning for microbiology and hematopathology slides
- TMA slide scanning
- Large format slide scanning
- True, 24-bit color data
- Standard TIFF file format

HIGH-PERFORMANCE VIEWING

- Instantly pan and zoom to any region and magnification, with thumbnail map of entire slide
- View multiple digital slides side-by-side
- Enhance images with real-time contrast, brightness, and gamma adjustments
- Annotate digital slides to provide descriptive information of regions of interest
- Create a path of views through a slide or group of slides by linking annotations and digital slides
- Extract selected regions of a digital slide to supplement presentations or publications

HEADQUARTERS

Aperio
1360 Park Center Drive
Vista, CA 92081

Main: 866.478.4111
Fax: 760.539.1116

info@aperio.com
support@aperio.com
www.aperio.com

EUROPEAN OPERATION

Aperio
3 The Sanctuary
Eden Office Park
Ham Green
Bristol BS20 0DD, UK

Main: 011 44 1275 375 123
Fax: 011 44 1275 373 501

europeinfo@aperio.com
europe-support@aperio.com
www.aperio.com

Aperio products are FDA cleared for specific clinical applications, and are intended for research use for other applications.

Copyright © 2008, Aperio. All rights reserved. ScanScope® and ImageScope™ are trademarks of Aperio.

IMC-1039-02072008

HOW TO GET STARTED:

STEP 1 Send us your glass slides. We will scan them using our award-winning ScanScope Scanner. Download the order form at Aperio.com, or contact us at 866.478.4111 to set an account.

STEP 2 Your scanned slides will be sent to you along with a DVD containing your digital slides. By request, our slide scanning service can include Internet access to your digital slides for 30 day free trial at slidehosting.com.

STEP 3 View, bookmark and annotate your digital slides. After a quick installation of our software, you will have complete functionality, affording you the advantage of powerful digital slide viewing, navigation, and analytic capabilities. If you choose to host your slides at www.slidehosting.com, you will be able to create bookmarks and annotations for your digital slides directly. We will provide you with a password to create and protect your files.

Call today and we'll scan up to 5 slides FREE, and we'll host your slides for up to 30 days FREE as well.

For information on how you can get started or to learn more about our solutions for digital pathology, including our family of ScanScope scanners and Spectrum™ Plus digital pathology management software visit www.aperio.com or contact our Inside Sales department at info@aperio.com or 866.478.4111.



bringing digital pathology to life